








SGE Inlet Liners with f.e.i.t.™

A revolutionary **NEW packing material technology** that has been designed to improve customers' experience in terms of liner longevity, analytical sensitivity, and liner-to-liner repeatability.

How long does f.e.i.t. last?

In a recent lifetime benchmarking study using the 8270 methodology with a soil matrix, **SGE inlet liner with f.e.i.t.** outperformed the competitive product by 4 times, achieving 400 injections before failure!

Why choose f.e.i.t.?

-  **Increased longevity**
-  **Improve data**
-  **Reduce cost and downtime**
-  **Resistant to dirty sample matrices**
-  **Touchless installation**

What does f.e.i.t. stand for:

f - folium l - liner
e - enhanced t - technology

Current portfolio offering:

f.e.i.t. is currently offered in the SGE deactivated bottom taper liner for GC instruments that utilize the industry standard geometry. Available in 1 and 5 pack configurations.



New biodegradable packaging to provide 100% recyclable material

**ALL YOU NEED IS
f.e.i.t.™**



SHOP NOW

Visit our website at [afplifsciences.com](https://www.afplifsciences.com)

Description	Pkg. 1	Pkg. 5
Taper Inlet Liner with f.e.i.t. with Pre-fitted o-ring, 4 mm ID, Length 78.5 mm	CL-100024	CL-100025