

Custom-made Packed Columns

Stainless Steel, Copper, Aluminum, Nickel and PTFE



Our custom-made OV[®]-Packed columns carry a 100% Satisfaction Guarantee!

Orders for custom-made columns will be given immediate attention and most will be shipped within 24 hours from receipt of order. We can supply packed columns made with any combination of support and liquid phase listed in this catalog. We will also make every effort to supply your requirements for other supports or liquid phases.

- Over 42 years experience in packing preparation and column techniques.
- 24 hour delivery.
- Complete with brass Swagelok fittings.
- All stainless steel tubing is our premium grade, manufactured for GC use only.
- Each column packed in individual container. (No lost ferrules or end caps)
- OV's guide to packings and columns included




Diagram 1




Diagram 2




Diagram 3




Diagram 4




Diagram 5




Diagram 6




Diagram 7





Diagram 8

Request A quote Example:



Request A Quote:

Stainless Steel, Copper, Aluminum, Nickel and PTFE Columns

Complete your information below and fax to (740) 373-9910. We will send you a firm quote.

	X							for	
↑		↑	↑	↑	↑	↑	↑		↑
Length		O. D. 1/8" 3/16" 1/4" 1/16"	Tubing Type Stainless Steel Copper Aluminum Nickel PTFE	Percent Loading	Liquid Phase	Solid Support	Mesh Size		GC Make & Model



SHOP NOW

Visit our website at
afplifsciences.com



Brass Fittings (standard)



Stainless Steel Fittings



Preconditioned

Company

Contact Names

Address

City

State

Zip

Phone

Fax

Email